

**Real-time-PCR Master Mix E5**  
**Cat-no: S185 - 100 rcs (2x1,25 ml)**

**Features:**

- The Master mix contains dUTP instead of dTTP
- The Master contains UNG (Uracil-N-Glycosylase)
- The Mix contains **ROX (100nM)** as passive Reference dye (it provides a baseline in multiplex reactions)
- It contains EvaGreen as fluorescent dye
- The qPCR / RTD-PCR Master mix E6 is ready-to-use and is optimized for high specificity and sensitivity because of optimized reaction buffer
- easy to us because ready-to-use Master Mix

**Applications:**

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Viral load determination

**Description:**

The Master Mix contains all reagents required for qPCR (except template and primer) in a premixed 2x concentrated ready-to-use solution. The high specificity and sensitivity of the mix is achieved by an optimized hot-start polymerase. Its activity is blocked at ambient temperature and switched on automatically at the onset of the initial denaturation. The thermal activation prevents the extension of non-specifically annealed primers and primer-dimer formations at low temperatures during PCR setup.

The mix offer dUTP instead of dTTP to prevent carry-over contaminations of DNA from previous PCR reactions.

**Concentration:** The Mastermix is 2x concentrated

**List of components qPCR / RTD-PCR Master mix:**

Hot-Start Polymerase for qPCR, dATP, dCTP, dGTP, dUTP, EvaGreen, ROX, UNG, optimized reaction buffer with KCl and MgCl<sub>2</sub>, stabilizers and enhancers,

**Transportation:** with blue ice

**Storage:** at 4°C for 3 months, at -20°C for more than 12 months, **Note:** protect from Light

**Usage:**

Components	Volume per reaction	final conc.
<b>2X qPCR / RTD-PCR Master mix E6, with UNG</b>	25 µl	1x
<b>Up-stream primer (10 µM stock)</b>	1,5 µl (range: 0,5-2.5 µl)	300 nM
<b>Down-stream primer (10µM stock)</b>	1,5 µl (range: 0.5-2,5 µl)	300 nM
<b>Template DNA</b>	5 µl (0.1-15 ng/ml plasmid DNA) (1-10 µg/ml genomic DNA)	< 500ng DNA
<b>Sterile dest. Water (included)</b>	up to 50 µl total reaction volume	

- vortex all solutions carefully before using and before PCR
- may you add the enzyme mix after Template DNA
- an individual optimization of annealing temperature may be necessary for new combinations of primers and Template DNA

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 25.10.2009 AS

Unless specified otherwise, all products of GeneON are sold for research use only.

*.. a good decision ..*

**General Thermo-Cycler protocol:**

**Note:** working with EvaGreen just select the optical setting for FAM or SYBR Green at the cycler

Step	Time	Temperature
UNG treatment	1x2 min	50°C
Initial denaturation	1-3 min	95°C
<b>30-40 Cycles:</b>		
Denaturation	15-30 sec	95°C
Annealing	30-65 sec	55-65°C
Extension	30 sec (per 500bp)	72-75°C

**Note:**

- an individual optimization of annealing temperature may be necessary for new combinations of primers and Template DNA

**Ordering information**

Cat.-no	Description	Amount
S185	qPCR / Real-time PCR Mastermix E6	100 rcs / 50µl PCR test vol.

*.. a good decision ..*