



Fluorescent DNA Loading dye

Maximo Fluorescent Reagent 6X Loading dye as substitute/replacement for Ethidium Bromide

Advantages for enduser:

- **Safety:** – Absence of mutagenicity and low toxicity (LC₅₀>5000mg/kg) as compared to Ethidium Bromide.
- **Low Environmental Impact:** Compliance with the Clean Water Act standards. No water pollution concern.
- **Sensitivity:** High degree of sensitivity as Ethidium Bromide.
- **Convenience:** Ready to Use; Same application **procedures as the 6X Loading Dye.**
- **Speed:** No de-staining requirement, low background value, and image displayed after coupling with the nucleic acid.
- **Compatibility:** Use the Blue Light or UV to detect the signal; Broad compatibility range.
- **Economic:** Non-hazardous product; No expenses required for the waste management.

Description:

Maximo Fluorescent Reagent is a non-mutagenic reagent that produces instant visualization of DNA bands upon Blue Light or UV illumination of agarose gels. Supplied in GeneON's 6X DNA Loading Buffer, Maximo Fluorescent Reagent is used to prepare DNA markers and samples for loading on agarose or polyacrylamide gels. Maximo Fluorescent Reagent is the most sensitive stain available for detecting the double-stranded DNA (dsDNA). It contains three tracking dyes (Bromophenol Blue, Xylene Cyanol FF, and Orange G) for visually tracking the DNA migration during the electrophoresis process. It is ideal for the environment requiring a safe, non-hazardous alternative to Ethidium Bromide.

Applications:

1. Vortex Maximo Fluorescent Reagent for 10 seconds prior to use.
 2. Dilute 1 part Maximo Fluorescent Reagent with 5 parts DNA sample and mix.
- Note: Maximo Fluorescent Reagent must be added to DNA markers in order to visualize the ladder bands simultaneously with the sample after electrophoresis.
3. Load sample and run according to standard procedures.
 4. After the electrophoresis, remove gel and place on UV or a visible-light transilluminator to immediately visualize bands.
 5. Gels can be post-stained with Ethidium Bromide if desired.

Tracking Dyes: Bromophenol Blue, Xylene Cyanol FF, and Orange G.

Storage and Transport:

Store at room temperature or at 4°C up to 12 months. For longer periods, store at -20°C.

Note: Maximo Fluorescent Reagent is light sensitive and should be stored protected from light.

Ordering information:

Cat.-no	Description	Amount
306-300	Maximo Fluorescent Reagent	1 ml
306-310	Maximo Fluorescent Reagent	5 x 1 ml

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..