



Phage Lambda DNA


Features:

λ DNA is used as one of the substrates in restriction enzymes research and for testing of restriction endonucleases activity.

Description/Preparation:

The double stranded DNA is isolated from bacteriophage lambda (cl857 ind1 Sam7). The molecular weight is 31.5×10^6 daltons and it is 48,502 base pairs in length.

The phage is isolated from the heat inducible lysogen E. coli λ cl857 S7 by gel filtration. The DNA is isolated from the purified phage by phenol/chloroform extraction and dialyzed against 10 mM Tris-HCl (pH 7.6) and 1 mM EDTA.



Concentration: 0.2-0.5 $\mu\text{g}/\mu\text{l}$

Storage Buffer: 10 mM Tris-HCl (pH 7.6), 1 mM EDTA.

Quality control: Gel analysis for purity. EcoRI and HindIII fragmentation patterns.

Transportation: Shipped on blue ice or room temperature

Storage: at -20°C for 24 months

Ordering information:

Cat.-no	Description	Amount
301-025	Phage Lambda DNA	500 μg
301-100	Phage Lambda DNA	5 x 500 μg

.. a good decision ..

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.taq-dna.com/>

Version: 25.10.2009 AS