



## MMLV Reverse Transcriptase

### Applications:

- RT PCR
- Synthesis of cDNA
- mRNA 5'-end Mapping by Primer Extension Analysis
- End-labeling of DNA
- Dideoxynucleotide Sequencing

### Description:

MMLV Reverse Transcriptase, encoded by Moloney Murine Leukemia Virus (MMLV RT) is an RNA-dependent DNA polymerase that synthesizes the cDNA first strand from a single-stranded RNA template to which a primer has been hybridized. MMLV RT will also extend primers hybridized to single-stranded DNA. Second strand cDNA synthesis can be achieved from some mRNA templates without an additional DNA polymerase.

**Concentration:** 200 u/μl

### Storage Buffer:

20 mM Tris-HCl, pH 7.4, 100 mM NaCl, 0,1mM EDTA, 50 % Glycerol, 0,1% IGEPAL, 1 mM DTT

### Reaction Buffer complete 10X: (Note: Composition and concentration has changed from 5x to 10x)

500 mM Tris-HCl (pH 8.3 at 25°C); 30 mM MgCl<sub>2</sub>; 750 mM KCl; 100 mM DTT.

**Dilution Buffer 1X:** 10 mM KH<sub>2</sub>PO<sub>4</sub> (pH 7.5); 0,1 mM EDTA; 200 mM NaCl; 7 mM 2-mercaptoethanol; 50% glycerol

### Unit definition:

One unit of the enzyme incorporates 1 nmol dTTP into acid-precipitable material in 10 minutes at 37°C, using poly(A) oligo dT as a template primer.

### Quality control:

**Endonuclease Activity:** 1 μg of Type 1 supercoiled plasmid DNA is incubated with 500 units of enzyme in 1X reaction buffer for one hour at 37°C. The supercoiled DNA is visualized on an ethidium bromide-stained agarose gel to verify absence of nicking or cutting.

**Nuclease Activity:** 50 ng of radio labelled DNA or RNA is incubated with 200 units of enzyme in 1X reaction buffer for one hour at 37°C, resulting in <1% release for both DNase and RNase.

**Purity:** >90% as judged by SDS-polyacrylamide gels with blue staining. MMLV RT is free of detectable RNase, and DNase (exo- and endonuclease) activities.

### Usage:

Standard Protocol:

We recommend to prepare 2 Mixes

### Mix I

Component	Amount/conc.
a. Total RNA or b. PolyA RNA	1-5 μg 50-500 ng
c. Strand-specific primer or d. oligo dT / random primer for each μg of RNA	10 pM 250-500 ng
sterile Water	up to 8 μl
<b>Incubation</b>	<b>Temperature</b>

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: [www.geneon.net](http://www.geneon.net) Version: 24-04-23 AS

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..



10 min	70 °C
10 - 15 min (for <b>c.</b> specific primers) or 5 min (for <b>d.</b> oligo dT / random primer )	room temperature  place on ice

### Mix II

Component	Amount/conc.
10X reaction buffer	2 µl
dNTP mix (10 mM of each = 40 mM)	1 µl
optional: RNasin	20-40 units
MMLV Reverse (200 u/µl)	200 units
sterile water	up to 20 µl
combine Mix I and Mix II and gently vortex	

Step	Temperature
30 - 115 min <sup>1.)</sup>	37 - 55°C <sup>2.)</sup>
10 min (Inactivation of enzyme)	65-70°C

1.) 30 min for cDNA with 500 bp; 115 min for 1,5 kb

2.) depends on the RNA: Higher temperatures (up to 55 °C) for higher structured RNA; Try to adjust the pH to 8.8

**Note:** MMLV used under standard 37°C reaction conditions is best for synthesis of 7 kb or less. Increasing reaction temp to 42°C may allow synthesis up to 10 kb. MMLV activity drops drastically above 42°C.

**Transportation:** on blue ice

**Storage:** at -20°C for 24 months

### Ordering information:

Cat.-no	Description	Amount
105-100	MMLV Reverse Transcriptase	10.000 units
105-250	MMLV Reverse Transcriptase	50.000 units

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

**E-Mail:** <mailto:info@geneon.net> **WEB:** [www.geneon.net](http://www.geneon.net) Version: 24-04-23 AS

Unless specified otherwise, all products of GeneON are sold for research use only.

*.. a good decision ..*