

Tween 20*

Nonionic detergent widely used in biochemical applications

Product	Tween 20
Quality	for molecular biology
Product number	8506
Formula	C ₅₈ H ₁₁₄ O ₂₆
Molecular weight	1227.72 g/mol
CAS number	9005-64-5
Storage	RT

*® Trademark of ICI America Inc.

Parameter	Specification
DNases/RNases/Proteases	not detectable
Density (d 20°C/4°C)	1.095 - 1.105
Hydroxyl-No.	approx. 100
Saponification No.	40 - 50
CMC (25°C)	5.9 x 10 ⁻⁵ mol/l
Classification	nonionic detergent
Solubility (25°)	100 g/L (H ₂ O)
Boiling point	> 100°C

Ordering information:

Cat.-no	Description	Amount
CE242	Tween 20	1000 ml

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..