



MATERIAL SAFETY DATA SHEET

**VIVANTIS TECHNOLOGIES SDN BHD
REVONGEN CORPORATION CENTER
12A, Jalan TP5, Taman Perindustrian UEP,
47600 Subang Jaya, Selangor, Malaysia.**

For Information:

Tel: +603-8025 1603

Fax: +603-80351637 / 1354

Email: info@vivantechnologies.com

URL: www.vivantechnologies.com

Document No.: MSDSrev03_GFPT

Date prepared: 1st August 2019

Reviewed: 6th August 2019

SECTION 1: CHEMICAL IDENTIFICATION

Code: COGF-FU-005; GF-FU-050; GF-FU-100

Name: GF-1 Fungus DNA Extraction Kit with 5 preps, 50 preps and 100 preps

Description: DNA extraction kit for a wide variety fungus tissues.

Intended Use: For research use and laboratory use only.

SECTION 2: COMPOSITION/INFORMATION ON INGREDIENTS

All concentrations are below the acceptable limits specific to each chemical.

Component(s)	Hazardous Ingredients (s)	CAS No	Concentration
FuL Buffer	a. Sodium Dodecyl Sulfate	151-21-3	≤1%
FuB Buffer	a. Guanidine Thiocyanate	593-84-0	≥30% <60%
Wash Buffer	No hazardous substance or mixture in the ingredients		
Elution Buffer	No hazardous substance or mixture in the ingredients		
Proteinase K	Proteinase, Tritirachium album serine	39450-01-6	>95%

Acute Effects

The ingredients in this diagnostic kit are irritants to skin, eyes and respiratory system. They are harmful if it comes into contact with the skin and toxic if ingested.

Chronic Effects

Hazardous in case of ingestion.

Pairing Nature with Scientific Discoveries

Vivantis Technologies Sdn Bhd 587399-D

A member of Revongen Corporation (SST No: B16-1808-21023214)

Revongen Corporation Center

No12A, Jalan TP5, Taman Perindustrian UEP,

47600 Subang Jaya, Selangor Darul Ehsan, Malaysia.

E : info@vivantechnologies.com

W: www.vivantechnologies.com

T: +6 03 8025 1603

F: +6 03 8025 1637

SECTION 3: HAZARDS IDENTIFICATION

GHS Classification

The extraction kit consists of individual ingredients for different components of kit.

Component(s)	Hazardous Ingredients (s)	CAS No	Concentration
FuL Buffer	a. Sodium Dodecyl Sulfate	151-21-3	≤1%
Classification (Regulation (EC) No1272/2008). Flammable Solid 2; H228 Acute Toxicity 4; H332 Skin Irritation 2; H315 Serious Eye Damage 1; H318 STOT – Respiratory System 3; H335 Chronic aquatic toxicity 3; H412			
Classification (67/548/EEC) No information available.			

Component(s)	Hazardous Ingredients (s)	CAS No	Concentration
FuB Buffer	a. Guanidine Thiocyanate	593-84-0	≥30% - <60%
Classification (Regulation (EC) No. 1272/2008) Acute Toxicity 4; H302 Acute Toxicity 4; H332 Hazardous to aquatic environment-chronic hazard 3; H412			
Classification (67/548/EEC) No information available.			

Component(s)	Hazardous Ingredients (s)
Wash Buffer	No hazardous substance or mixture in the ingredients
Classification (Regulation (EC) No. 1272/2008) Not a hazardous substance or mixture.	
Classification (67/548/EEC) Not a hazardous substance or mixture.	

Component(s)	Hazardous Ingredients (s)
Elution Buffer	No hazardous substance or mixture in the ingredients
Classification (Regulation (EC) No. 1272/2008) Not a hazardous substance or mixture.	
Classification (67/548/EEC) Not a hazardous substance or mixture.	

Component(s)	Hazardous Ingredients (s)	CAS No	Concentration
Proteinase K	Proteinase, Triticirachium album serine	39450-01-6	>95%
Classification (Regulation (EC) No. 1272/2008) Skin Irritation 2; H315 Eye Irritation 2; H319 Respiratory Sensitivity 1; H334 Skin Sensitivity 1; H317 STOT Single Exposure 3; H335			
Classification (67/548/EEC) Xi; R36/37/38/42 Xi; Xi; R36/37/38 Xn; R42 Xi; R43			

Classification of the whole kit according to CLASS regulation 2013

Acute toxicity, Category 4, Oral; H302: Harmful if swallowed
 Acute toxicity, Category 2, Inhalation; H332: Harmful if inhaled
 Acute toxicity, Category 2, Dermal; H317: May cause an allergic skin reaction
 STOT Single Exposure, Respiratory tract irritation, Category 3; H335: May cause respiratory irritation
 Chronic aquatic toxicity, Category 3; H412: Harmful to aquatic life with long lasting effects

Label elements

Labeling according to regulation (EC) No 1272/2008

Hazard pictograms



Pairing Nature with Scientific Discoveries

Signal word

Danger

Hazard statements

H228	Flammable solid
H302 + H312 + H332	Harmful if swallowed or if inhaled
H314	Causes severe skin burns and eye damage
H315	Causes skin irritation
H317	May cause an allergic skin reaction
H318	Causes serious eye damage
H319	Causes serious eye irritation
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled
H335	May cause respiratory irritation
H412	Harmful to aquatic life with long lasting effects

Precautionary statements

Prevention:

P210	Keep away from heat
P261	Avoid breathing dust / fume / gas / mist / vapours / spray
P264	Wash skin thoroughly after handling
P270	Do not eat, drink or smoke when using this product
P273	Avoid release to the environment
P280	Wear eye protection / face protection
P280	Wear protective gloves
P284	Wear respiratory protection

Response:

P301 + P330 + P331	If swallowed, rinse mouth. Do not induce vomiting.
P302 + P352	If on skin, wash with plenty of soap and water.
P304 + P340 + P312	If inhaled, remove person to fresh air and keep comfortable for breathing. Call a poison center or doctor/physician if you feel unwell.
P305 + P351 + P338 + P310	If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center or doctor / physician.
P342 + P311	If experiencing respiratory symptoms, call a poison center or doctor / physician.
P362 + P364	Take off contaminated clothing and wash it before reuse.
P501	Dispose of contents/container to an approved waste disposal plant.

Other hazards

Not applicable

SECTION 4: FIRST-AID MEASURES

In case of contact with eyes, immediately flush with copious amounts of water for at least 15 minutes. Remove contact lenses, if present and easy to do.

In case of contact with skin, immediately wash with soap and copious amounts of water for several minutes. Take off immediately all contaminated clothing.

In case of ingestion, wash out mouth with water provided the person is conscious. Do not induce vomiting. Never give anything by mouth to an unconscious person. Call a physician immediately. Do not attempt to neutralize.

In case of inhalation, move to fresh air. Give artificial respiration if not breathing. Call a physician.

Notes to physician are that the first aid procedure should be established in consultation with the doctor responsible for industrial medicine.

SECTION 5: FIRE FIGHTING MEASURES

Extinguishing Media

Suitable: Water spray, carbon dioxide, dry chemical powder or alcohol-resistant foam.

Unsuitable: High volume water jet.

Do not allow run-off from fire fighting to enter drains or water courses.

Development of hazardous combustion gases or vapours possible in the event of fire. Fire may cause evolution of sulphur oxides and nitrogen oxides.

Wear self contained breathing apparatus and protective clothing and equipment to prevent contact with skin and eyes for firefighting.

Collect contaminated fire extinguishing water separately. This must not be discharged into drains. Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations.

SECTION 6: ACCIDENTAL RELEASE MEASURES

Personal Precautions

Wear the protective clothing, rubber gloves and mask. Wear self-contained breathing apparatus, rubber boots and goggles if necessary.

Do not breathe vapours, aerosols. Avoid substance contact. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Environmental Precautions

Do not let product enter drains. Prevent further leakage or spillage if safe to do so.

If the product contaminates rivers and lakes or drains inform respective authorities.

Measures for Cleaning / Collecting

Soak up with inert absorbent material: sand, silica gel, acid binder, universal binder, sawdust.

Collect spilled liquid with liquid-binding material or inert absorbent and place in closed container for disposal.

Dispose of in accordance with federal, state and local environmental regulations.

Wash spill site after material pickup is complete.

SECTION 7: HANDLING AND STORAGE

Handling

Good laboratory techniques should be used when handling. Do not eat, drink, smoke or apply cosmetics in laboratory areas. Do not pipette reagents or samples by mouth. Use reagents according to the product insert.

Work under hood or open-air space. Do not inhale substance / mixture. Avoid generation of vapours / aerosols.

Apply preventive skin protection. Change the contaminated clothing immediately after work. Wash hands and face after working with substance.

Avoid extreme temperatures during transport.

Storage

Store all components as directed in the product insert. Keep tightly closed in a dry and well ventilated place.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

Wear appropriate personal protective equipment when working with components or samples, including laboratory coats, disposable gloves and eye protection. Avoid hand and mouth contact. Wash hands as soon as possible after handling components or samples.

Handle in accordance with good industrial hygiene and safety practice.

Control Parameters

Exposure Limits: Contains no substances with occupational exposure limit values
Engineering measures: Ensure adequate ventilation, especially in confined areas

Predicted No Effect Concentration (PNEC) according to Regulation (EC) No. 1907/2006

Sodium Dodecyl Sulphate

Environmental Compartment	Value
Fresh Water	0.176 mg/L
Marine Water	0.0176 mg/L
Fresh Water Sediment	6.97 mg/L
Marine Sediment	0.697 mg/L
Soil	1.29 mg/L

Exposure Controls

Respiratory Protection

In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand Protection

Wear suitable protective gloves. Gloves material: compatible chemical-resistant gloves; nitrile rubber gloves.

Pairing Nature with Scientific Discoveries

Eye & Face Protection	Eye wash bottle with pure water. Tightly fitting safety goggles. Wear face-shield and protective suit for abnormal processing problems.
Skin & Body Protection	Wear suitable protective clothing and impervious clothing. Choose body protection according to the amount and concentration of the dangerous substance at the work place.

Environmental Exposure Controls

Do not let product enter drains.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Information for PL Buffer, PB Buffer, Wash Buffer, Elution Buffer, Proteinase K

Appearance:	Liquid, Colorless
Odor:	No odor
pH:	FuL Buffer: pH8-9.5 FuB Buffer: pH5.5 Wash Buffer: pH7-9 Elution Buffer: pH8-9 Proteinase K: pH7.5
Melting Point:	No information available
Boiling Point:	No information available
Flash Point:	No information available
Autoignition Temperature:	No information available
Decomposition Temperature:	No information available
Evaporation rate:	No information available
Flammability:	No information available
Upper Explosion Limit:	No information available
Lower Explosion Limit:	No information available
Vapour Pressure:	No information available
Relative Density:	No information available
Specific Gravity:	No information available
Water Solubility:	All buffers soluble
Partition coefficient: n-octanol/water:	No information available
Explosive Properties:	Not classified as explosive
Oxidizing Properties:	Not classified as oxidizing
Other information:	No information available

Pairing Nature with Scientific Discoveries

Vivantis Technologies Sdn Bhd 587399-D

A member of Revongen Corporation (SST No: B16-1808-21023214)

Revongen Corporation Center

No12A, Jalan TP5, Taman Perindustrian UEP,
47600 Subang Jaya, Selangor Darul Ehsan, Malaysia.

E : info@vivantechnologies.com
W : www.vivantechnologies.com

T : +6 03 8025 1603
F : +6 03 8025 1637

SECTION 10: STABILITY AND REACTIVITY

The contents of the kit are stable under normal handling and storage conditions as stated in the product insert until the expiration date indicated on the corresponding label.

Reactivity: No dangerous reaction known under conditions of normal use
 Possibility of hazardous reactions: Reacts with oxidizing agents; keep away from heat and sources of ignition. Risk of dust explosion
 No decomposition if stored and applied as directed

Conditions to avoid: Extremes of temperature and direct sunlight. Heat, flames and sparks

Incompatible materials: No dangerous reaction known under conditions of normal use
 Materials to avoid: strong acids, oxidizing agents, bases

Polymerization: Hazardous polymerization does not occur

Hazardous decomposition products:

In case of fire hazardous decomposition products may be produced such as: carbon monoxide, carbon dioxide, and unburned hydrocarbons (smoke), nitrogen oxides, ammonia, gaseous hydrogen chloride, carbon oxides, sulphur oxides and sodium oxides.

SECTION 11: TOXICOLOGICAL INFORMATION

Information for Sodium Dodecyl Sulfate (FuL Buffer)

Acute oral toxicity: LD50 Oral (Rat): 1.200 mg/kg
 Acute inhalation toxicity: No information available
 Acute dermal toxicity: LD50 Dermal (Rat, male and female): >2.001 mg/kg

Skin corrosion / irritation: Causes skin irritation
 Serious eye damage / eye irritation: Causes serious eye damage
 Respiratory or skin sensitization: Not classified based on available information
 Germ cell mutagenicity: Not classified based on available information
 Carcinogenicity effects: Not classified based on available information
 Reproductive toxicity: Not classified based on available information
 Aspiration toxicity: Not classified based on available information
 STOT – Single exposure: May cause respiratory irritation.
 STOT – Repeated exposure: Not classified based on available information

Information for Guanidine Thiocyanate (FuB Buffer)

Acute Toxicity

Acute oral toxicity: LD50 Oral (Rat): 593 mg/kg
 Acute inhalation toxicity: Acute toxicity estimate: 1.5 mg/L; Test atmosphere: dust/mist
 Acute dermal toxicity: Acute toxicity estimate: 1100 mg/kg

Pairing Nature with Scientific Discoveries

Skin corrosion / irritation:	No data available
Serious eye damage / eye irritation:	No data available
Respiratory or skin sensitization:	No data available
Germ cell mutagenicity:	No data available
Carcinogenicity effects:	No data available
Reproductive toxicity:	Not classified based on available information
Aspiration toxicity:	Not classified based on available information
STOT – Single exposure:	Not classified based on available information
STOT – Repeated exposure:	Not classified based on available information

Information for Wash Buffer and Elution Buffer

Acute Toxicity:	Not classified based on available information
Skin corrosion / irritation:	Not classified based on available information
Serious eye damage / eye irritation:	Not classified based on available information
Respiratory or skin sensitization:	Not classified based on available information
Germ cell mutagenicity:	Not classified based on available information
Carcinogenicity effects:	Not classified based on available information
Reproductive toxicity:	Not classified based on available information
Aspiration toxicity:	Not classified based on available information
STOT – Single exposure:	Not classified based on available information
STOT – Repeated exposure:	Not classified based on available information

Information for Proteinase K

Acute Toxicity:	Not classified based on available information
Skin corrosion / irritation:	Causes skin irritation and / or dermatitis
Serious eye damage / eye irritation:	Causes eye irritation; may cause irreversible eye damage
Respiratory or skin sensitization:	May cause skin allergy or may cause sensitization by inhalation
Germ cell mutagenicity:	Not classified based on available information
Carcinogenicity effects:	Not classified based on available information
Reproductive toxicity:	Not classified based on available information
Aspiration toxicity:	Not classified based on available information
STOT – Single exposure:	May cause respiratory irritation.
STOT – Repeated exposure:	Not classified as specific target organ toxicant, repeated exposure.

Additional toxicological information

The product is not subject to classification according to the calculation method of the General EU Classification Guidelines for preparation as issued in the latest version. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

SECTION 12: ECOLOGICAL INFORMATION

The extraction kit contains **Sodium Dodecyl Sulfate (PL Buffer)** and **Guanidine Thiocyanate (PB Buffer)** known to be hazardous to the environment or not degradable in waste water treatment plants. Self assessment: Hazardous for water. Do not allow product or large quantities of it to reach ground water, water course or sewage system.

Toxicity

Sodium Dodecyl Sulfate

Fish: LC50 (Leuciscus idus (Golden orfe)): 29 mg/L

LC50 (Lepomis macrochirus (Bluegill sunfish)): 4.1 mg/L; Exposure time: 96hr

LC50 (Oncorhynchus mykiss (rainbow trout)): 3.6 mg/L; Exposure time: 96hr

NOEC (Pimephales promelas (fathead minnow)): 1.36 mg/L; Exposure time: 42d

Daphnia and other aquatic invertebrates:

EC50 (Daphnia magna (Water flea)): 5.55 mg/L; Exposure time: 48hr

NOEC (Ceriodaphnia dubia (water flea)): 0.88 mg/L; Exposure time: 7d

Microorganisms:

LC50 (activated sludge): 135 mg/L; Exposure time: 3hr

Chronic aquatic toxicity: Harmful to aquatic life with long lasting effects.

Toxicity Data on Soil: Not expected to adsorb on soil.

Other organisms relevant to the environment: No data available.

Guanidine Thiocyanate

Fish: LC50 (Poecilia reticula (guppy)): 89.1 mg/L; Exposure time: 96 hr

Daphnia and other aquatic invertebrates: EC50 (Daphia (water flea)): 42.4 mg/L; Exposure time: 48hr

Chronic aquatic toxicity:

Harmful to aquatic life with long lasting effects

Toxicity data on soil:

Not expected to absorb on soil

Other organisms relevant to the environment:

No data available

Persistence and Degradability

Sodium Dodecyl Sulfate

Biodegradation: 95%; 28d;

Bioaccumulative Potential

Sodium Dodecyl Sulfate

Bioaccumulation: Cyprinus carpio (Carp); 3d; 3.9-5.3 (Bioconcentration Factor)

Partition coefficient n-octanol/water: log Pow: -2.03 (20°C)

Guanidine Thiocyanate

Partition coefficient n-octanol/water: log Pow: -1.38

Mobility in Soil

No data available for the extraction kit, including Sodium Dodecyl Sulfate and Guanidine Thiocyanate.

Results of PBT and vPvB assessment

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Other adverse effects

No information available.

SECTION 13: DISPOSAL CONSIDERATIONS

Do not contaminate ponds, waterways or drains with chemical or used container. Used components from this product should be disposed to a licensed waste management company.

Waste material must be disposed of in accordance with federal, state and local environmental regulations.

SECTION 14: TRANSPORT INFORMATION

IATA / ADR / RID / ADN / DOT-US / IMDG

Not regulated as dangerous product in the meaning of transport regulations.

UN Number:	Not regulated as dangerous product
UN proper shipping name:	Not regulated as dangerous product
Transport hazard class:	Not regulated as dangerous product
Packing group:	Not regulated as dangerous product
Environmental Hazards:	Not regulated as dangerous product
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code:	Not regulated as dangerous product

SECTION 15: REGULATORY INFORMATION

International Inventories

REACH:	Not comply
CH INV:	Complies
TSCA:	Complies
DSL:	Complies
AICS:	Complies
NZIoC:	Complies
ENC:	Not comply
ISHL:	Not comply
KECI:	Complies
PICCS:	Complies
IECES:	Complies



- R 22: Harmful if swallowed.
R 36/37/38: Irritating to eyes, respiratory system and skin.
R 36/37: Irritating to eyes and skin.
R 41: Risk of serious damage to eyes.
R42: May cause sensitization by inhalation.
R43: May cause sensitization by skin contact.
S 26: Rinse eyes immediately with plenty of water.
S 36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

Korea / Taiwan / China

No information available

US Federal Regulations

SARA 313

Not regulated

US State Regulations

California Proposition 65

Not regulated

Canada Regulatory Information

WHMIS Classification

This product has been classified in accordance with the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

DSL: Yes

NDSL: No

SECTION 16: OTHER INFORMATION

The information contained in this MSDS relates only to the material(s) designed and does not relate to use(s) in combination with any other material, process(es) and /or chemical reaction(s). Vivantis Technologies Sdn. Bhd. provides this information in good faith, from sources believed to be accurate; however, Vivantis assumes no liability for its accuracy or completeness, and thus shall not be held liable for any damage resulting from handling or from contact with the above product.

All Vivantis products are supplied for manufacturing, research and laboratory use only. Researchers and laboratory personnel intending to use any of these products for medical investigation on human are solely responsible for such use and for compliance with the pertinent regulations of the United States Food & Drug Administration (US-FDA) and other regulations. We do not assume liability for damages resulting from the handling, use and/or disposal of these products, from their use in violation of patent or other rights or reliance upon this information.

Pairing Nature with Scientific Discoveries

Vivantis Technologies Sdn Bhd 587399-D

A member of Revongen Corporation (SST No: B16-1808-21023214)

Revongen Corporation Center

No12A, Jalan TP5, Taman Perindustrian UEP,
47600 Subang Jaya, Selangor Darul Ehsan, Malaysia.

E : info@vivantechnologies.com
W: www.vivantechnologies.com

T: +6 03 8025 1603
F: +6 03 8025 1637